



HIGH DEFINITION OPTICS



HIGH DEFINITION OPTICS



Continuing our drive of optical optimization ConMed Linvatec offers the highest quality optics throughout the whole imaging chain. HD Laparoscopes offer improved optical resolution at the perimeter of the visual field and increased depth of field while maintaining an exceptional level of illumination. The HD lens design ensures maximized transmission of high contrast light from the laparoscope to the camera head.

- Improved light transmission
- Improved resolution around the perimeter of the visual field
- Improved depth of field



10mm HD Autoclavable Laparoscope	Eyecup	Cartridge
0 Degree	HD1000	HD1001
30 Degree	HD1030	HD1031
45 Degree	HD1044	HD1045

5.5mm HD Autoclavable Laparoscope		
0 Degree	HD5500	HD5501
30 Degree	HD5530	HD5531
45 Degree	HD5544	HD5545

5.5mm Extended Length HD Autoclavable Laparoscope		
0 Degree	HD5500XL	HD5501XL
30 Degree	HD5530XL	HD5531XL
45 Degree	HD5544XL	HD5545XL

Laparoscope Specifications:

	HD 10mm		HD 5.5mm		HD 5.5mm XL
Outside Diameter	10mm	Outside Diameter	5.5mm	Outside Diameter	5.5mm
Angle of View	0, 30, 45	Angle of View	0, 30, 45	Angle of View	0, 30, 45
Field of View	75	Field of View	75	Field of View	75
Working Length	328mm	Working Length	300mm	Working Length	425mm



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com